# scil v-Parvo

## Immunological Rapid Test


**FOR VETERINARY USE ONLY!**

### TEST INFORMATION

The scil v-Parvo is a diagnostic test for animals (dogs and cats) and is used for the detection of specific antigens of Parvovirus (CPV, FPV) in faeces. Canine Parvovirus type 2 (CPV2) is the causative agent of gastroenteritis. The new types of CPV2 are also infectious for felines and can be transmitted between the species. A negative result with the scil v-Parvo antigen test does not completely rule out a parvovirus infection: Faecal shedding of parvovirus antigen at detectable levels typically occurs between 3-12 days post-exposure and usually correlates with the onset of clinical signs. In puppies with moderate maternal antibody levels, viral shedding may be delayed by 1-2 days relatively to the onset of clinical signs. Virus shedding begins to wane by day 8-10 (post infection). The scil v-Parvo is an antigen detecting test that confirms suspected cases of parvovirus infection within 10 minutes with the possibility to start the treatment directly on site.

### TEST COMPONENTS

**Note:** Prior to the test usage, the reaction field shows a green line in the control line region. This is a quality indicator and will be washed away by the sample fluid during the test procedure.

### PLEASE NOTE PRIOR TO USE

Please use a new test cartridge for every individual test as cartridges are for single-use only. scil Rapid Test kits are for veterinary use only. Use only test components provided by scil animal care company. Use the test cassette within 60 minutes after opening the pouch and place the test cassette in a horizontal position on a smooth surface while the test is performed. Note the amount of sample material needed. A too high amount of feces may disturb the test procedure. After opening the pouch, use the test cassette within one hour. Consider the test results as invalid after the read out time. Do not use the test after the expiration date on the pouch. Fecal samples may be infectious. Dispose all contaminated materials properly and disinfect the work area after the test execution.

### STORAGE


scil Rapid Test kit should be stored between 2-30°C.

### SAMPLE MATERIAL

Best sample material are fresh dog or cat faeces. Pooled faecal samples are possible. Then the sample must be stirred well prior to testing. If it is impossible to test on the day of sample collection, the faecal sample can be stored max. 6 days at 2-8°C. The faecal sample must have obtained room temperature (18-25°C) before test execution. **Stir the faecal sample thoroughly prior to testing.** Avoid picking up cat litter or other coarse particles.


### CORRECT SAMPLING

Use the swab to stir the faecal sample well. Wipe off excessive faecal material on the walls of the faecal collection tube, so that the tip of the swab is coated with faecal material.

### TEST PROCEDURE

Open the sample tube. Take the cotton swab and put the cotton swab with the faecal sample into the sample tube, containing the reagent. Stir up well the fluid with the cotton swab so that the sample material dissolves. The reagent should become considerably discolored.

Close the sample tube with the reagent tightly. The reagent in the sample tube will treat and conserve the sample. Shake the sample tube well. Open the aluminum pouch of the test cassette and place the cassette on a horizontal surface. Break up the tip of the sample tube by pressing it firmly. Hold the sample tube with the lid down over the sample well of the test cassette.


Apply three (3) drops to the sample well of the test cassette. If the liquid does not run up the test strip after a few seconds, add another drop to the sample well until the liquid starts running.

**Note:** If too much faeces or mucus got into the sample well, so that the liquid does not run up the test strip after a few seconds, you can press with the upper tip of the cotton swab into the sample well to reactivate the run of the test.

### TEST EVALUATION

The result of the test can be read after 10 minutes.

For a positive result, two red lines appear in the reaction field of the test cassette. **A red line in the T-region (T) of the reaction field indicates a positive test result (parvovirus antigen detected).** A faint test line is also considered as positive result. The second red line in the C-region (C) indicates the control line, which is consistent with the correct performance of the test. The C-line is not a reference line and may have a different line intensity than the T-Line.

### Invalid Result:

If no control line appears after 10 minutes, the test is invalid. In this case, the test was maybe not properly conducted or the test is expired. If this occurs, a new test must be conducted.

An excessive amount of sample faeces can cause a brown line in the T-region (T), which must not be seen as a positive result. The test is invalid and should be repeated.

### TEST PERFORMANCE

	Sensitivity	Specificity	Reference	n
CPV/ FPV	93.33%	99.99%	ELISA*	100

\*Enzyme-linked Immunosorbent Assay

### REFERENCE

**Decaro N. et al.:** 2005, Clinical and virological findings in pups naturally infected by canine parvovirus type 2 Glu-426 mutant. J Vet Diagn Invest 17:133-138.  
**Green Tree, Ernst G., Schimke, Ernst:** 2007, Clinic of the dog diseases, Enke publishing house, pp. 111 ff.  
**Greene CE, ed.:** 2006, Infectious diseases of the dog and cat. Elsevier Saunders, Oxford, UK.  
**Schmitz S., et al.:** 2009, Brief research report. J Vet Diagn Invest 21: 344-345


### MANUFACTURER

produced in Germany for scil animal care company GmbH.  
 Dina-Weissmann-Allee 6, 68519 Viernheim, Germany

Tel.: +49 (0) 6204 78 90 0,  
 E-Mail: [info-de@scilvet.com](mailto:info-de@scilvet.com)

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Prodotto in Germania per scil animal care company GmbH.